## Quick Reference Summary Sheet TECH 03a - Intravenous Tail Vein Injections in Adult Mice

## **Procedure:**

- 1. Pre-warm sterile substance to be administered. Do not exceed 37°C.
- 2. Weigh animal. Calculate volume to administer.
  - a. Bolus Injection (over 1-2 seconds) = 5ml/kg
  - b. Slow Bolus Injection (over 5-6 seconds) = 10ml/kg
- 3. Fill syringe and needle with pre-warmed substance. Maintain sterility. Ensure no air bubbles.
  - a. Mice: 0.3-1.0 ml syringe, 27 30 G needle
- 4. Safely warm animal. Do not leave unattended. Do not overheat.
- 5. Place animal in restrainer or anesthetize.
- 6. Identify lateral tail vein on either side of tail (Figure 1; grey arrows).
- Occlude vessel at base of tail (optional). Apply gentle traction to keep tail straight.
  - a. Swab tail with isopropyl alcohol-moistened gauze
  - b. Gently bend tail over ring finger (Figure 2)
- 8. Insert needle, parallel to vein, in middle third of tail (Figure 3)
  - a. Pull back on plunger to look for blood in needle hub (optional)
  - b. Release occlusion of vein (if doing)
- 9. Inject substance and look for the vein to clear of blood.
  - a. If unsuccessful (blanched skin and bleb; Figure 4), another attempt can be made closer to the mouse or the tail vein on the other side. Two attempts per tail vein may be performed.
- 10. Leave needle in tail for 2-3 seconds so the substance can be cleared from the vein.
- 11. Remove needle. Immediately apply pressure to the insertion site for ~10 seconds until bleeding has stopped.
- 12. Discard syringe and needle into sharps container. Do not recap.
- 13. Remove animal from restrainer or recover from anesthesia and examine again for any bleeding.
- 14. Return animal to its cage and observe for any complications.
- 15. Note procedure on cage card/monitoring records.

## **Potential Complications:**

- · Local reaction:
  - E.g. bleeding, bruising, burns, necrosis or swelling
- Systemic reaction:
  - E.g. overheating, shock, dyspnea, weight loss, hunched, piloerection, lethargy, isolated, pain







